

Oak Medical Services Ltd.

TQs Electronic Tourniquet

Product Specifications:

- Main LCD: 14cm x 14cm, displays 2 channels: Cuff Pressure, Pre-set Pressure Timer. Alarms signs.
- Secondary LCD: User settings i.e. Surtrak timer, Pressure / Procedure Logs & Engineering Menu and readouts.
- Cuff inflation from 0 – 760MMhg Paediatric to XL Cuff quick inflation.
- Digital display with high grade membrane touch button controls - Not Touch Screen as this is not favoured in the UK market because of accidental activation.
- TQs has 2 channels that can be used separately (perform bilateral / Biers Block procedures).
- Automatic timer and Surtrak programable surgical time tracking with alarm.
- Battery back-up is 4 hours in use from full charge with healthy battery.
- Internal Battery charger.
- Mains Cable 110 - 240V
- Audio and visual alarms: Faults, Cuff Leaks,
- Colour coded armoured cuff supply hoses.
- Facia Cuff Hose holders to help prevent broken cuff connectors.
- Rear PC interface for engineer investigation and readout.
- Rear Earth point for Testing
- Easy to read user guide / instruction /diagram on unit.